



Safety data sheet according to UK REACH

Printing date 14.03.2025

Version number 1

Revision: 14.03.2025

1 Identification of the substance/mixture and of the company/undertaking

· **Product identifier**

· **Trade name:** DermaDam™

· **Article number:** SDS 181-001.06R01, 311, 314, S311

· **CAS Number:**

9006-04-6

· **EC number:**

232-689-0

· **Relevant identified uses of the substance or mixture and uses advised against**

No further relevant information available.

· **Application of the substance / the mixture**

The dental dam is a single use, non-invasive, non-sterile device that is designed to be used only by qualified dental professionals for the described purpose. It is not a measuring device. It is designed for use by the professional in a medical environment for periods of up to 30 minutes.

· **Uses advised against**

It is the responsibility of the dental practitioner to determine if the use of dental dam is appropriate for each individual patient. The dental dam should not be used in situations where the skin in the relevant area is broken, the patient has known sensitivity or other clinical considerations indicate that it should not be used.

· **Details of the supplier of the safety data sheet**

· **Manufacturer/Supplier:**

Manufactured by Four D Rubber Company Ltd for:

Ultradent Products Inc.

505 W. Ultradent Drive (10200 S)

South Jordan, UT 84095-3942

USA

onlineordersupport@ultradent.com

EC Responsible Person

Ultradent Products GmbH

Am Westhover Berg 30

51149 Cologne Germany

Email: infoDE@ultradent.com

Office Phone: +49(0)2203-35-92-0

· **Further information obtainable from:** Customer Service

· **Emergency telephone number:**

CHEMTREC (NORTH AMERICA) : +1 (800) 424-9300

(INTERNATIONAL) : +(703) 527-3887

2 Hazards identification

· **Classification of the substance or mixture**

· **Classification according to Regulation (EC) No 1272/2008**

The substance is not classified, according to the GB CLP regulation.

· **Label elements**

· **Labelling according to Regulation (EC) No 1272/2008** Void

· **Hazard pictograms** Void

· **Signal word** Void

(Contd. on page 2)

Safety data sheet according to UK REACH

Printing date 14.03.2025

Version number 1

Revision: 14.03.2025

Trade name: **DermaDam™**

(Contd. of page 1)

· **Hazard statements** *Void*

3 Composition/information on ingredients

- **Substances**
- **Description:** *Mixture of substances listed below with nonhazardous additions.*
9006-04-6 Natural rubber latex
- **Identification number(s)**
- **EC number:** 232-689-0

4 First aid measures

- **Description of first aid measures**
- **General information:** *No special measures required.*
- **After inhalation:**
Not applicable under normal temperature conditions. Remove from exposure in the event of inhalation of fumes/vapours at elevated temperatures and seek medical advice.
- **After skin contact:**
Not considered hazardous for short term exposure. Product may cause type IV allergic reactions on contact within susceptible individuals. Remove exposure, perform first aid if necessary and seek medical advice in the event of type IV allergic reaction to the product.
- **After eye contact:** *If product contacts eye, rinse with sterile eye wash.*
- **After swallowing:** *If product is ingested, perform first-aid measures, and seek medical advice.*
- **Most important symptoms and effects, both acute and delayed** *No further relevant information available.*
- **Indication of any immediate medical attention and special treatment needed**
No further relevant information available.

5 Firefighting measures

- **Extinguishing media**
- **Suitable extinguishing agents:**
Water
Foam
Carbon dioxide
Use fire extinguishing methods suitable to surrounding conditions.
- **Special hazards arising from the substance or mixture** *No further relevant information available.*
- **Advice for firefighters:**
- **Protective equipment:** *No special measures required.*
- **Additional information**
Will ignite if in contact with a naked flame. Observe non smoking regulations in storage areas.

6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures** *Not required.*
- **Environmental precautions:** *Do not allow to enter sewers/ surface or ground water.*
- **Methods and material for containment and cleaning up:** *Pick up mechanically.*
- **Reference to other sections**
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.

(Contd. on page 3)

Safety data sheet according to UK REACH

Printing date 14.03.2025

Version number 1

Revision: 14.03.2025

Trade name: **DermaDam™**

(Contd. of page 2)

See Section 13 for disposal information.

7 Handling and storage

- **Precautions for safe handling:**
See product labeling.
No special measures required.
- **Information about fire - and explosion protection:** No special measures required.
- **Conditions for safe storage, including any incompatibilities**
Ideally store below 26C and away from heaters. If stored below 15C, rubber may become brittle and lose elasticity. Avoid damp storage conditions, store below 70% relative humidity. Avoid exposure to ionising radiation sources. Store away from sources of ozone such as mercury vapour lamps and high voltage electrical equipment. Avoid contact with copper, manganese and associated alloys containing either metal which will cause staining. Avoid contact with oils, solvents and greases, which will destroy synthetic rubber.
- **Storage:**
- **Requirements to be met by storerooms and receptacles:**
Product should be stored flat in supplied container with lid closed when not in use. Protect from prolonged exposure to light, especially bright sunlight or strong U.V. light.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** None.
- **Specific end use(s)** No further relevant information available.

8 Exposure controls/personal protection

- **Control parameters**
- **Ingredients with limit values that require monitoring at the workplace:** Not required.
- **Additional information:** The lists valid during the making were used as basis.
- **Exposure controls**
- **Appropriate engineering controls** No further data; see section 7.
- **Individual protection measures, such as personal protective equipment**
- **General protective and hygienic measures:**
The usual precautionary measures are to be adhered to when handling chemicals.
- **Respiratory protection:** Not required.
- **Hand protection**
The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.
Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.
Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation
- **Material of gloves**
The product is intended to be used with latex or non-latex gloves by the operator. The product should not be used on broken skin.
The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer.
- **Penetration time of glove material**
The exact breakthrough time has to be found out by the manufacturer of the protective gloves and has to be observed.
- **Eye/face protection** Not required.
- **Body protection:** Protective work clothing

GB

(Contd. on page 4)

Safety data sheet

according to UK REACH

Printing date 14.03.2025

Version number 1

Revision: 14.03.2025

Trade name: **DermaDam™**

(Contd. of page 3)

9 Physical and chemical properties

· Information on basic physical and chemical properties

· General Information

· Physical state	Solid
· Colour:	Not determined.
· Odour:	Characteristic 'rubber' odour
· Odour threshold:	Not determined.
· Melting point/freezing point:	Undetermined.
· Boiling point or initial boiling point and boiling range	Undetermined.
· Flammability	Product is not flammable.
· Lower and upper explosion limit	
· Lower:	Not determined.
· Upper:	Not determined.
· Flash point:	Not applicable.
· Decomposition temperature:	Not determined.
· pH	Not applicable.
· Viscosity:	
· Kinematic viscosity	Not applicable.
· Dynamic:	Not applicable.
· Solubility	
· water:	Insoluble.
· Partition coefficient n-octanol/water (log value)	Not determined.
· Vapour pressure:	Not applicable.
· Density and/or relative density	
· Density at 20 °C:	0.00321 g/cm ³
· Relative density	Not determined.
· Vapour density	Not applicable.

· Other information

· Appearance:	
· Form:	Rubber Sheeting, either as a roll or cut squares
· Important information on protection of health and environment, and on safety.	
· Ignition temperature:	Not determined.
· Explosive properties:	Product does not present an explosion hazard.
· Change in condition	
· Evaporation rate	Not applicable.

· Information with regard to physical hazard classes

· Explosives	Void
· Flammable gases	Void
· Aerosols	Void
· Oxidising gases	Void
· Gases under pressure	Void
· Flammable liquids	Void
· Flammable solids	Void
· Self-reactive substances and mixtures	Void
· Pyrophoric liquids	Void
· Pyrophoric solids	Void
· Self-heating substances and mixtures	Void
· Substances and mixtures, which emit flammable gases in contact with water	Void
· Oxidising liquids	Void
· Oxidising solids	Void

(Contd. on page 5)

Safety data sheet

according to UK REACH

Printing date 14.03.2025

Version number 1

Revision: 14.03.2025

Trade name: **DermaDam™**

(Contd. of page 4)

· Organic peroxides	Void
· Corrosive to metals	Void
· Desensitised explosives	Void

10 Stability and reactivity

- **Reactivity** No further relevant information available.
- **Chemical stability**
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **Possibility of hazardous reactions:** No dangerous reactions known.
- **Conditions to avoid:** No further relevant information available.
- **Incompatible materials:** No further relevant information available.
- **Hazardous decomposition products:** No dangerous decomposition products known.

11 Toxicological information

- **Information on hazard classes as defined in Regulation (EC) No 1272/2008**
- **Acute toxicity** Based on available data, the classification criteria are not met.
- **Primary irritant effect:**
- **Skin corrosion/irritation**
Not thought to cause skin irritation. Traces of the processing additives within synthetic rubber latex can cause contact irritant dermatitis in some cases, leading to redness, soreness, dryness or cracking of the skin.
Based on available data, the classification criteria are not met.
- **Serious eye damage/irritation** Based on available data, the classification criteria are not met.
- **Respiratory or skin sensitisation** Based on available data, the classification criteria are not met.
- **Germ cell mutagenicity** Based on available data, the classification criteria are not met.
- **Carcinogenicity** Based on available data, the classification criteria are not met.
- **Reproductive toxicity** Based on available data, the classification criteria are not met.
- **STOT-single exposure** Based on available data, the classification criteria are not met.
- **STOT-repeated exposure** Based on available data, the classification criteria are not met.
- **Aspiration hazard**
Product is intended to be placed in the mouth, improper use could cause an aspiration hazard. Product is intended to be used by professionals, and the patient should be continually supervised when the product is in use.
Based on available data, the classification criteria are not met.
- **Additional toxicological information:**
Allergies: Traces of the processing additives within the synthetic rubber latex can give rise to type IV allergic reactions, leading to redness, soreness, dryness or cracking of the skin.
- **Information on other hazards**

- **Endocrine disrupting properties**

Substance is not listed.

12 Ecological information

- **Toxicity**
- **Aquatic toxicity:** No further relevant information available.
- **Persistence and degradability** biodegradable
- **Other information:** The product is biodegradable.
- **Bioaccumulative potential** No further relevant information available.
- **Mobility in soil** No further relevant information available.

(Contd. on page 6)

Safety data sheet according to UK REACH

Printing date 14.03.2025

Version number 1

Revision: 14.03.2025

Trade name: **DermaDam™**

(Contd. of page 5)

- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.
- **Endocrine disrupting properties** The product does not contain substances with endocrine disrupting properties.
- **Other adverse effects**
- **Additional ecological information:**
- **General notes:**
*Water hazard class 2 (German Regulation) (Self-assessment): hazardous for water
 Do not allow product to reach ground water, water course or sewage system.
 Danger to drinking water if even small quantities leak into the ground.*

13 Disposal considerations

- **Waste treatment methods**
- **Recommendation**
Dispose of contents/container in accordance with international, federal, state, and local regulations.
- **Uncleaned packaging:**
- **Recommendation:** *Disposal must be made according to official regulations.*

14 Transport information

<ul style="list-style-type: none"> · UN number or ID number · ADR, IMDG, IATA 	<i>not regulated</i>
<ul style="list-style-type: none"> · UN proper shipping name · ADR, ADN, IMDG, IATA 	<i>not regulated</i>
<ul style="list-style-type: none"> · Transport hazard class(es) · ADR, ADN, IMDG, IATA · Class 	<i>not regulated</i>
<ul style="list-style-type: none"> · Packing group · ADR, IMDG, IATA 	<i>not regulated</i>
<ul style="list-style-type: none"> · Environmental hazards: 	<i>Not applicable.</i>
<ul style="list-style-type: none"> · Special precautions for user 	<i>Not Applicable</i>
<ul style="list-style-type: none"> · Maritime transport in bulk according to IMO instruments 	<i>Not applicable.</i>
<ul style="list-style-type: none"> · UN "Model Regulation": 	<i>not regulated</i>

15 Regulatory information

<ul style="list-style-type: none"> · Safety, health and environmental regulations/legislation specific for the substance or mixture 	
<ul style="list-style-type: none"> · NIOSH-Ca (National Institute for Occupational Safety and Health) 	<i>Substance is not listed.</i>
<ul style="list-style-type: none"> · Poisons Act 	
<ul style="list-style-type: none"> · Regulated explosives precursors 	<i>Substance is not listed.</i>

(Contd. on page 7)

Safety data sheet according to UK REACH

Printing date 14.03.2025

Version number 1

Revision: 14.03.2025

Trade name: DermaDam™

(Contd. of page 6)

· **Regulated poisons**

Substance is not listed.

· **Reportable explosives precursors**

Substance is not listed.

· **Reportable poisons**

Substance is not listed.

· **Directive 2012/18/EU**

· **Named dangerous substances - ANNEX I** Substance is not listed.

· **Chemical safety assessment:** A chemical safety assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· **Department issuing SDS:** Environmental, Health, and Safety

· **Contact:** Customer Service

· **Abbreviations and acronyms:**

ADR: Accord relatif au transport international des marchandises dangereuses par route (European Agreement Concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonised System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

NIOSH: National Institute for Occupational Safety

GB