

Printing date 03.03.2025 Version number 1 Revision: 03.03.2025

1 Identification of the substance/mixture and of the company/undertaking

- · Product identifier
- · Trade name: UltraTempTM REZ II Catalyst (Regular, Fast)
- · Article number: SDS 439-001.03R01, 1007965, 1007604
- Relevant identified uses of the substance or mixture and uses advised against Professional dental luting & filling agent
- · Application of the substance / the mixture Professional dental luting & filling agent
- · Details of the supplier of the safety data sheet
- · Manufacturer/Supplier:

Ultradent Products Inc.

505 W. Ultradent Drive (10200 S)

South Jordan, UT 84095-3942

USA

onlineordersupport@ultradent.com

EC Responsible Person

Ultradent Products GmbH

Am Westhover Berg 30

51149 Cologne Germany

Email: infoDE@ultradent.com

Emergency Phone: +49(0)2203-35-92-0

Ultradent Products Inc.

505 W. Ultradent Drive (10200 S)

South Jordan, UT 84095-3942

USA

onlineordersupport@ultradent.com

- · Further information obtainable from: Customer Service
- · Emergency telephone number:

CHEMTREC (NORTH AMERICA) : (800) 424-9300 (INTERNATIONAL) : +(703) 527-3887

2 Hazards identification

- · Classification of the substance or mixture
- · Classification according to Regulation (EC) No 1272/2008



environment

Aquatic Acute 1 H400 Very toxic to aquatic life.

Aquatic Chronic 1 H410 Very toxic to aquatic life with long lasting effects.



Skin Sens. 1

H317 May cause an allergic skin reaction.

- · Label elements
- · Labelling according to Regulation (EC) No 1272/2008 Void
- · Hazard pictograms GHS07, GHS09

(Contd. on page 2)

Printing date 03.03.2025 Version number 1 Revision: 03.03.2025

Trade name: UltraTempTM REZ II - Catalyst (Regular, Fast)

(Contd. of page 1)

- · Signal word Warning
- Hazard-determining components of labelling:

Triethylene Glycol Dimethacrylate

· Hazard statements

H317 May cause an allergic skin reaction.

H410 Very toxic to aquatic life with long lasting effects.

· Precautionary statements

P101 If medical advice is needed, have product container or label at hand.

P102 Keep out of reach of children.

P103 Read carefully and follow all instructions.

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

P273 Avoid release to the environment.

P280 Wear protective gloves.

P333+P313 If skin irritation or rash occurs: Get medical advice/attention.

P321 Specific treatment (see on this label).

P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

3 Composition/information on ingredients

- · Mixtures
- Description: Mixture of substances listed below with nonhazardous additions.

· Dangerous components:					
CAS: 1314-13-2	Zinc Oxide	≥25-<40%			
EINECS: 215-222-5	🕸 Aquatic Acute 1, H400; Aquatic Chronic 1, H410				
CAS: 109-16-0	Triethylene Glycol Dimethacrylate	>10-≤25%			
EINECS: 203-652-6	♦ Skin Sens. 1, H317				
	Dimethyl-p-toluidine	≥0.1-<5%			
EINECS: 202-805-4	Acute Tox. 3, H301; Acute Tox. 3, H311; Acute Tox. 3, H331; STOT RE 2, H373; Aquatic Chronic 3, H412				

[·] Additional information: For the wording of the listed hazard phrases refer to section 16.

4 First aid measures

- · Description of first aid measures
- · General information: Immediately remove any clothing soiled by the product.
- · After inhalation:

Supply fresh air and to be sure call for a doctor.

In case of unconsciousness place patient stably in side position for transportation.

- · After skin contact: Immediately wash with water and soap and rinse thoroughly.
- · After eye contact: Rinse opened eye for several minutes under running water.
- · After swallowing: If symptoms persist consult doctor.
- · Most important symptoms and effects, both acute and delayed No further relevant information available.
- · Indication of any immediate medical attention and special treatment needed

No further relevant information available.

5 Firefighting measures

- · Extinguishing media
- · Suitable extinguishing agents: Use fire extinguishing methods suitable to surrounding conditions.
- · Special hazards arising from the substance or mixture No further relevant information available.

(Contd. on page 3)

Printing date 03.03.2025 Version number 1 Revision: 03.03.2025

Trade name: UltraTempTM REZ II - Catalyst (Regular, Fast)

(Contd. of page 2)

· Advice for firefighters:

· Protective equipment: No special measures required.

6 Accidental release measures

- · Personal precautions, protective equipment and emergency procedures Not required.
- · Environmental precautions:

Do not allow product to reach sewage system or any water course.

Inform respective authorities in case of seepage into water course or sewage system.

Do not allow to enter sewers/surface or ground water.

· Methods and material for containment and cleaning up:

Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).

Dispose contaminated material as waste according to section 13.

Ensure adequate ventilation.

· Reference to other sections

See Section 7 for information on safe handling.

See Section 8 for information on personal protection equipment.

See Section 13 for disposal information.

7 Handling and storage

· Precautions for safe handling:

Ensure good ventilation/exhaustion at the workplace.

Prevent formation of aerosols.

- · Information about fire and explosion protection: No special measures required.
- · Conditions for safe storage, including any incompatibilities
- · Storage:
- · Requirements to be met by storerooms and receptacles: No special requirements.
- · Information about storage in one common storage facility: Not required.
- · Further information about storage conditions: See product labelling.
- · Specific end use(s) Professional dental luting & filling agent

8 Exposure controls/personal protection

- · Control parameters
- · Ingredients with limit values that require monitoring at the workplace:

The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.

- · Additional information: The lists valid during the making were used as basis.
- · Exposure controls
- · Appropriate engineering controls No further data; see section 7.
- · Individual protection measures, such as personal protective equipment
- · General protective and hygienic measures:

Immediately remove all soiled and contaminated clothing.

Wash hands before breaks and at the end of work.

· Respiratory protection:

In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use self-contained respiratory protective device.

(Contd. on page 4)

Printing date 03.03.2025 Version number 1 Revision: 03.03.2025

Trade name: UltraTempTM REZ II - Catalyst (Regular, Fast)

(Contd. of page 3)

· Hand protection



The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.

Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

· Material of gloves

The selection of suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

· Penetration time of glove material

The exact breakthrough time has to be found out by the manufacturer of the protective gloves and has to be observed.

- · Eye/face protection Goggles recommended during refilling
- · Body protection: Protective work clothing

9 Physical and chemical properties

· Information on basic physical and chemical properties

· General Information

· Physical state Fluid

· Colour: According to product specification

· Odour: Characteristic · Odour threshold: Not determined. · Melting point/freezing point: Undetermined. · Boiling point or initial boiling point and boiling range >250 °C

· Flammability Not applicable.

· Lower and upper explosion limit

· Lower: Not determined. · Upper: Not determined. · Flash point: >150 °C · Auto-ignition temperature: 255 °C

· Decomposition temperature: *Not determined.*

Not applicable (non-aqueous)

 $\cdot pH$

· Viscosity:

· Kinematic viscosity Not determined. · Dynamic: Not determined.

·Solubility

Not miscible or difficult to mix.

· Partition coefficient n-octanol/water (log value) Not determined.

· Vapour pressure at 20 °C: 0 hPa

· Density and/or relative density

· Density at 20 °C: 1.65-1.75 g/cm3 · Relative density Not determined. · Vapour density Not determined.

(Contd. on page 5)

Printing date 03.03.2025 Version number 1 Revision: 03.03.2025

Trade name: UltraTempTM REZ II - Catalyst (Regular, Fast)

(Contd. of page 4)

٠	Otl	ier	in	form	ation

· Appearance:

Form: Viscous liquid

Important information on protection of health and environment, and on safety.

• Ignition temperature: Product is not selfigniting.

• Explosive properties: Product does not present an explosion hazard.

· Change in condition

• Evaporation rate Not determined.

Information with regard to physical hazard classes

Void · Explosives Void · Flammable gases · Aerosols Void · Oxidising gases Void · Gases under pressure Void · Flammable liquids Void · Flammable solids Void · Self-reactive substances and mixtures Void · Pyrophoric liquids Void · Pyrophoric solids Void · Self-heating substances and mixtures Void · Substances and mixtures, which emit flammable gases Void in contact with water Oxidising liquids Void · Oxidising solids Void · Organic peroxides Void Void · Corrosive to metals

10 Stability and reactivity

· Desensitised explosives

- · Reactivity No further relevant information available.
- · Chemical stability
- Thermal decomposition / conditions to be avoided: No decomposition if used according to specifications.

Void

- · **Possibility of hazardous reactions:** No dangerous reactions known.
- · Conditions to avoid: No further relevant information available.
- · Incompatible materials: No further relevant information available.
- · Hazardous decomposition products: No dangerous decomposition products known.

11 Toxicological information

- Information on hazard classes as defined in Regulation (EC) No 1272/2008
- · Acute toxicity Based on available data, the classification criteria are not met.

LD/LC50 values relevant for classification:

ATE (Acute Toxicity Estimates)

Oral	LD50	100,000 mg/kg
Inhalative	LC50/4 h	1,400 mg/l

(Contd. on page 6)

Printing date 03.03.2025 Revision: 03.03.2025 Version number 1

Trade name: UltraTempTM REZ II - Catalyst (Regular, Fast)

	(Contd. of page
Zinc Oxide	
LD50	7,950 mg/kg (mouse)
	>5,000 mg/kg (rat)
LC50 Fish	>320 mg/l (Fish)
LC50/4 h	>5.7 mg/l (mouse)
LC50(Daphnia magna) (static)	0.098 mg/l (daphnia)
riethylene Glycol Dimethacryla	nte
LD50	>5,000 mg/kg (rat)
LC50 Fish	16.4 mg/l (Fish) (Toxicity to fish)
LD50	>2,000 mg/kg (mouse)
methyl-p-toluidine	
LC50 Fish	46-52 mg/l (Fish)
LD50	>2,000 mg/kg (rabbit)
LC50/4 h	1.4 mg/l (rat)
	LC50 Fish LC50/4 h LC50/4 h LC50(Daphnia magna) (static) Triethylene Glycol Dimethacryla LD50 LC50 Fish LD50 methyl-p-toluidine LC50 Fish LD50

- · Information on other hazards
- · Endocrine disrupting properties None of the ingredients is listed.

12 Ecological information

· Toxicity

· Aquatic toxicity:			
1314-13-2 Zinc Oxide			
EC50	2.066 mg/kg (Fish)		
109-16-0 Triethylene Glycol Dimethacrylate			
EC50	>100 mg/kg (Algae)		
Biodegradability 28 days (Aerobic) (Biodegradability testing)			
Aqua toxicity	32 mg/l (daphnia) (No Observed Effect Concentration)		

- · Persistence and degradability No further relevant information available.
- · Bioaccumulative potential No further relevant information available.
- · Mobility in soil No further relevant information available.
- · Results of PBT and vPvB assessment
- · **PBT:** Not applicable.
- · vPvB: Not applicable.
- · Endocrine disrupting properties The product does not contain substances with endocrine disrupting properties.
- · Other adverse effects
- · Remark: Very toxic for fish
- · Additional ecological information:
- · General notes:

Water hazard class 2 (German Regulation) (Self-assessment): hazardous for water

Do not allow product to reach ground water, water course or sewage system.

Danger to drinking water if even small quantities leak into the ground.

Also poisonous for fish and plankton in water bodies.

(Contd. on page 7)

Printing date 03.03.2025 Version number 1 Revision: 03.03.2025

Trade name: UltraTempTM REZ II - Catalyst (Regular, Fast)

Very toxic for aquatic organisms

(Contd. of page 6)

13 Disposal considerations

- · Waste treatment methods
- · Recommendation

Dispose of contents/container in accordance with international, federal, state, and local regulations.

- · Uncleaned packaging:
- · Recommendation: Disposal must be made according to official regulations.

1 1 7	7	 C .	. •
	ransport	OVMAI	701

· UN number or ID number		UN	number	or ID	number
--------------------------	--	----	--------	-------	--------

- · ADR, IMDG, IATA not regulated
- · UN proper shipping name
- · ADR, ADN, IMDG, IATA not regulated
- · Transport hazard class(es)
- · ADR, ADN, IMDG, IATA
- · Class not regulated
- · Packing group
- · ADR, IMDG, IATA not regulated
- Environmental hazards: Not applicable.
- · Special precautions for user Not Applicable
- · Maritime transport in bulk according to IMO

instruments Not applicable.

· UN "Model Regulation": not regulated

15 Regulatory information

- Safety, health and environmental regulations/legislation specific for the substance or mixture
- NIOSH-Ca (National Institute for Occupational Safety and Health)

None of the ingredients is listed.

- · Poisons Act
- · Regulated explosives precursors

None of the ingredients is listed.

· Regulated poisons

None of the ingredients is listed.

· Reportable explosives precursors

None of the ingredients is listed.

· Reportable poisons

None of the ingredients is listed.

(Contd. on page 8)

Printing date 03.03.2025 Version number 1 Revision: 03.03.2025

Trade name: UltraTempTM REZ II - Catalyst (Regular, Fast)

(Contd. of page 7)

- · Directive 2012/18/EU
- · Named dangerous substances ANNEX I None of the ingredients is listed.
- · Seveso category E1 Hazardous to the Aquatic Environment
- Qualifying quantity (tonnes) for the application of lower-tier requirements 100 t
- Qualifying quantity (tonnes) for the application of upper-tier requirements 200 t
- · Chemical safety assessment:

Device is biocompatible when used as directed by dental professionals per ISO 10993-1

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· Relevant phrases from Section 3

H301 Toxic if swallowed.

H311 Toxic in contact with skin.

H317 May cause an allergic skin reaction.

H331 Toxic if inhaled.

H373 May cause damage to organs through prolonged or repeated exposure.

H400 Very toxic to aquatic life.

H410 Very toxic to aquatic life with long lasting effects.

H412 Harmful to aquatic life with long lasting effects.

- · Department issuing SDS: Environmental, Health, and Safety
- · Contact: Customer Service
- · Abbreviations and acronyms:

ADR: Accord relatif au transport international des marchandises dangereuses par route (European Agreement Concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonised System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

NIOSH: National Institute for Occupational Safety

ATE: Acute toxicity estimate values

Acute Tox. 3: Acute toxicity - Category 3

Skin Sens. 1: Skin sensitisation – Category 1

STOT RE 2: Specific target organ toxicity (repeated exposure) – Category 2

Aquatic Acute 1: Hazardous to the aquatic environment - acute aquatic hazard - Category 1

Aquatic Chronic 1: Hazardous to the aquatic environment - long-term aquatic hazard - Category 1

Aquatic Chronic 3: Hazardous to the aquatic environment - long-term aquatic hazard – Category 3

* * Data compared to the previous version altered.